

AVISCERA BIOSCIENCE

Human Asprosin / Fibrillin-1 (2732-2871) Recombinant Biotinylated

01-50B Description

A DNA sequence encoding the human Asprosin /Fibrillin-1 (Ser²⁷³²-His²⁸⁷¹) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was biotinylated with water soluble Biotin and dialysis against PBS, pH 7.4.

Formulation

Lyophilized 50 μ g human Asprosin/ Fibrillin-1 (2732-2871) Biotinylated in 100 μ l of PBS, Carry free.

Reconstitution & Storage

Add 100 μ l PBS to the vial to prepare a working stock solution at 500 μ g/mL. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at - 20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: human Asprosin / Fibrillin-1 (Ser²⁷³²-His²⁸⁷¹)

Product Information

Code 00229-01-50B

Name Human Asprosin Biotinylated

Lot No.

Size 50 μg

Species Human

Sequence Ser2732-2871His

Protein ID NP_000129.3

Gene ID NM 000138.4

MW 19.7 KD

Tag His tag on N terminal

Source E. Coli

Purity >95% in SDS-PAGE gel

PBS lyophilized Form

without preservatives

Carry Free

Storage $-20 \,^{\circ}$ C $^{\sim}$ -70 $^{\circ}$ C

Reconstitution 100 μ l

AVISCERA BIOSCIENCE

2348 Walsh Ave, Suite C Santa Clara, CA 95051

USA

Tel: (408) 982 0300 Fax: (408) 982 0301

Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.