

AVISCERA BIOSCIENCE

Anti-Human Endothelial Cell-Specific Molecule 1 (ESM1) Monoclonal Antibody

Product Information

Code A00318-05-100

Name Human ESM1

Mab

Clone No. 95-1B2

Lot No.

Size 100 µl

Species Human

Host Mice

Immunogen ESM1 (Human),

rec.

Ab Type IgG

Purification sequential

precipitation

Formulation lyophilized

Form without preservatives

preservativ

Carry free

Storage -20 ° C

Specificity Human only

Reconstitutio

PBS, 100 μl

Application ELISA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, *E. coli*-derived, recombinant human ESM1. That antibody was purified by sequential precipitation caprylic acid and ammonium sulphate.

Formulation

100 μg of Anti Human ESM1 IgG in 100 μl of PBS without preservatives was lyophilized.

Reconstitution

Add 100 μ l of PBS to the vial to prepare antibody stock solution at 100 μ l. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 μ L per vial) and stored frozen at -20° C to -70° C in a manual defrost freezer for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. Avoid repeated freeze-thaw cycles.

Specificity

This antibody has been selected for its ability to recognize human ESM1 on direct ELISAs.

Applications

Indirect ELISA - This antibody can be used at 3.34ng/mL with the appropriate secondary reagents to detect human ESM1 on indirect ELISA.

Optimal dilutions should be determined by each laboratory for each application.

ORDER INFORMATION AVISCERA BIOSCIENCE, INC

2348 WALSH AVE., SUITE C SANTA CLARA, CA 95051 TEL: +001 408 982 0300

Info@Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.