

## **AVISCERA BIOSCIENCE**

# Rabbit Anti- Anterior Gradient 2 (AGR2) (Human) IgG

Code A00529—01-

100

Name Anti AGR2 (H)

IgG

Clone No. N/A

**Product Information** 

Lot No.

Size 100 μg

Species Human

Host Rabbit

Immunogen AGR2 (Human) mature form

Ab Type IgG

Purification Protein A

Lyophilized

Formulation Form without

Preservatives

Carry Free

Storage -20 ° C

Specificity Human

Application IHC, WB

**ELISA** 

AVISCERA BIOSCIENCE

2348 Walsh Ave., Suite C Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Info@AvisceraBioscience.com www.AvisceraBioscience.com

### Preparation

This antibody was produced from a rabbit immunized with purified recombinant human Anterior Gradient 2 (AGR2). That IgG was purified by Protein A affinity.

#### **Formulation**

100  $\mu g$  of purified Anti Anterior Gradient 2 (AGR2) (Human) IgG in 100  $\mu l$  of PBS without preservatives was lyophilized.

#### Reconstitution

Add 100  $\mu$ l of PBS to the vial to prepare antibody stock solution at 100  $\mu$ g/100  $\mu$ l. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted ( by 10  $\mu$ L per vial) and stored frozen at -20° C to -70° C in a manual defrost freezer for up six months without detectable loss of activity.

#### Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.** 

#### Specificity

This antibody has been selected for its ability to recognize human AGR2 on indirect ELISAs, immunohistochemistry and western blot.

#### **Applications**

Indirect ELISA - This antibody can be used at 0.02  $^{\sim}$ 0.04  $\mu g$  /mL with the appropriate secondary reagents to detect recombinant human AGR2 on indirect ELISA.

Immunohistochemistry-That Antibody can be used at 2  $^{\sim}4~\mu g$  /mL with the appropriate secondary antibody to detect AGR2 in paraffin embedded human lung , breast and pancreas cancer tissues slides (ABC).

Western Blot-This antibody can be used at 1-2  $\mu$ g /mL with the appropriate secondary antibody in western blot analysis to detect recombinant human AGR2.

Optimal dilutions should be determined by each laboratory for each application.

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.