



Anti Human Vaspin Monoclonal Antibody

Product Information

Code	A00560-02-100
Name	Anti Vaspin (Human) MAB
Clone No.	1H12E11
Lot No.	20110972
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human Vaspin rec.
Ab Type	IgG
Purification	Protein G Lyophilized
Formulation	Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

Preparation

This antibody was produced from an hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, E. coli-derived, recombinant human Vaspin. That antibody was purified by Protein G affinity.

Formulation

100 µg of Anti Human Vaspin Monoclonal Antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (1000µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for 12 months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize recombinant human Vaspin in indirect ELISAs.

Applications

Indirect ELISA Assay - This antibody can be used at 0.0625 ~0.125 µg/ml to detect recombinant human vaspin on indirect ELISA.

Optimal dilutions should be determined by each laboratory for each application.

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.