



AVISCERA BIOSCIENCE

Anti Human CTRP5 (Complement C1q Tumor necrosis factor-Related Protein 5) Monoclonal Antibody

Product Information

Code	A00594-08-100
Name	Human CTRP5 Mab
Clone No.	F9
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human CTRP5, Full length rec.
Ab Type	Purified IgG
Purification	Protein G affinity
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified recombinant human CTRP5, Full Length, His Tag on N-Terminal. That antibody was purified Protein G affinity.

Formulation

100 µg of Anti human CTRP5 Monoclonal antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (1000ug/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human CTRP5 Full Length on indirect ELISAs.

Applications

Indirect ELISA - This antibody can be used at 0.031-0.062 µg/ml to detect human CTRP5 Full Length on indirectly ELISA.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.