



AVISCERA BIOSCIENCE

ANTI-HUMAN BETA DEFENSIN 101 IGG

Product Information

Code	A00858-06-100
Name	Anti human BD101 IgG
Clone No.	
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	Human BD101 full length
Ab Type	IgG
Purification	Protein A
Formulation	lyophilized
Carry	N/A
Storage	-20 ° C
Specificity	Human
Application	ELISA

Preparation

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant human BD101 full length. That antibody was purified with Protein A affinity.

Formulation

100 µg of rabbit anti-human BD101 IgG in 100 µl of PBS without preservatives was lyophilized.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human BD101 full length and mature form on western blot and elisa.

Applications

Western Blot - This antibody can be used as antibody at 1-5 µg/ml with the appropriate detection antibody to detect human BD101 on transfected human tumor cells .

Direct ELISA - This antibody can be used at 1: 4000 (0.25 µg/mL) with the appropriate secondary reagents to detect recombinant human BD101 full length and mature form on indirectly elisa.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE, INC

2348 Walsh Ave, Suite C,
Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Fax: (408) 982 0301

info@Aviscerabioscience.com

www.Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.